



## APPLICATION

- Cavity preparation
- Crown & Bridge preparation
- Final modification of composite restorations
- Trimming and contouring on the dental materials: ceramics, precious metal, nonprecious metals, model cast alloys acrylics, etc.

## SAFETY RECOMMENDATIONS

1. Be sure that the bur is inserted fully and tightly into the handpiece.
2. Select suitable speed according to recommended speed table.
3. Never exceed maximum speed.
4. The recommended work pressure is 0.3 - 2N (30-200p) for FG type.  
The recommended work pressure is 2 - 5 N (200-500p) for HP type.
5. Avoid high contact pressure and leverage.
6. Dentists: please always use spray cooling (minimum 50ml/min)
7. The Recommended Speed Table (FG type) :

Round




<b>100</b>	<b>101</b>	<b>102</b>	<b>103</b>	<b>104</b>	<b>105</b>	<b>106</b>	<b>360</b>	<b>361</b>	<b>362</b>	<b>410</b>	<b>105L</b>
801/008	801/009	801/010	801/012	801/014	801/016	801/018	801/021	801/023	801/025	801/031	801L/016
001/008	001/009	001/010	001/012	001/014	001/016	001/018	001/021	001/023	001/025	001/031	697/016



Inverted with Collar




<b>131</b>	<b>132</b>	<b>133</b>	<b>134</b>	<b>135</b>	<b>136</b>
806/009	806/010	806/012	806/014	806/016	806/018
019/009	019/010	019/012	019/014	019/016	019/018



Round with Collar



<b>111</b>	<b>112</b>	<b>113</b>	<b>114</b>	<b>115</b>	<b>116</b>	<b>381</b>
802/009	802/010	802/012	802/014	802/016	802/018	802/023
002/009	002/010	002/012	002/014	002/016	002/018	002/023



Inverted Cone



<b>120</b>	<b>121</b>	<b>122</b>	<b>123</b>	<b>124</b>	<b>125</b>	<b>126</b>	<b>384</b>
805/008	805/009	805/010	805/012	805/014	805/016	805/018	805/022
010/008	010/009	010/010	010/012	010/014	010/016	010/018	010/022



Ogival Cylinder



<b>216</b>	<b>217</b>	<b>218</b>	<b>315</b>	<b>316</b>	<b>317</b>	<b>318A</b>	<b>318</b>	<b>319</b>
876/010	877/010	877/012	878/010	878/012	878/014	879/012	879/014	879/016
287/010	288/010	288/012	289/010	289/012	289/014	290/012	290/014	290/016



Conical Cylinder



<b>256</b>	<b>257</b>	<b>352A</b>	<b>352</b>	<b>352B</b>	<b>357A</b>	<b>357</b>
884/010	884/012	885/010	885/012	885/014	886/012	886/014
129/010	129/012	130/010	130/012	130/014	131/012	131/014



Cylinder



<b>201</b>	<b>201B</b>	<b>202</b>	<b>203</b>	<b>204</b>	<b>205</b>	<b>206</b>	<b>206B</b>	<b>206C</b>	<b>206D</b>	<b>301</b>	<b>302</b>	<b>302B</b>	<b>303</b>	<b>304</b>	<b>304B</b>
835/008	835/009	835/010	835L/010	835/012	835/014	836/012	836/014	836/016	836/018	837/012	837/014	837/016	837/018	842/012	842/014
108/008	108/009	108/010	108L/010	108/012	108/014	110/012	110/014	110/016	110/018	111/012	111/014	111/016	111/018	112/012	112/014

Round






<b>304C</b>	<b>276</b>	<b>431</b>	<b>303A</b>	<b>433</b>
842/016	836KR/018	836KR/023	837KR/018	837KR/023
112/016	157/018	157/023	158/018	158/023

Round End Cylinder



<b>270</b>	<b>270B</b>	<b>271</b>	<b>272</b>	<b>273</b>	<b>347</b>	<b>348</b>	<b>349</b>	<b>349B</b>	<b>349C</b>	<b>350</b>	<b>351</b>
838/008	838/009	838/010	838/010	838/012	880/010	880/012	880/014	880/016	880/018	881/012	881/014
138/008	138/009	138/010	139/010	139/012	140/010	140/012	140/014	140/016	140/018	141/012	141/014





**Incanted Cone**

<b>404</b>	<b>417</b>	<b>324</b>	<b>324B</b>	<b>325</b>	<b>326</b>	<b>326B</b>	<b>406</b>
856/021	856/025	850/012	850/014	850/016	850/018	850B/018	850/021
198/021	198/025	199/012	199/014	199/016	199/018	199B/018	199/021

•   •   •   •   •   •   •   •

**Ogival End Taper (Gingival Curettage)**

<b>220</b>	<b>221</b>	<b>222</b>	<b>223</b>	<b>370</b>	<b>371</b>	<b>372</b>	<b>372B</b>	<b>372C</b>	<b>405</b>	<b>373A</b>
877K/012	877K/014	877K/016	877K/018	878K/010	878K/012	878K/014	878K/016	878K/018	878K/021	879K/010
297/012	297/014	297/016	297/018	298/010	298/012	298/014	298/016	298/018	298/021	299/010

•   •   •   •   •   •   •   •   •   •   •

**Flat End Taper (Truncated Cone)**

<b>214</b>	<b>407</b>	<b>305</b>	<b>306</b>	<b>307</b>	<b>308</b>	<b>401</b>	<b>409</b>	<b>310</b>	<b>311</b>	<b>312</b>	<b>313</b>	<b>402</b>	<b>411</b>	<b>505</b>	<b>403</b>	<b>263</b>	<b>463</b>	<b>355</b>
846L/018	846L/025	847/012	847/014	847/016	847/018	847/021	847/025	848/012	848/014	848/016	848/018	848/021	848/025	848/031	848L/018	846KR/018	846KR/023	847KR/018
171L/018	171L/025	172/012	172/014	172/016	172/018	172/021	172/025	173/012	173/014	173/016	173/018	173/021	173/025	173/031	174/018	545/018	845/023	546/018

**Round End Taper (Domed Tr)**

<b>455</b>	<b>243</b>	<b>244</b>	<b>245</b>	<b>246</b>	<b>250</b>	<b>251</b>	<b>252</b>	<b>253</b>	<b>320A</b>	<b>320</b>	<b>321</b>	<b>322</b>	<b>323</b>	<b>323B</b>
847KR/023	849/009	849/010	849/012	849/014	855/010	855/012	855/014	855/016	856A/010	856/012	856/014	856/016	856/018	856B/018
546/023	196/009	196/010	196/012	196/014	197/010	197/012	197/014	197/016	198A/010	198/012	198/014	198/016	198/018	198B/018

•   •   •   •   •   •   •   •   •   •   •   •   •   •

**Diabolo (Amalgam Remover)**

<b>140</b>	<b>141</b>	<b>142</b>	<b>143</b>
813/010	813/012	813/016	813/018
032/010	032/012	032/016	032/018

**Inverted Cone Long**

<b>226</b>	<b>227</b>	<b>228</b>	<b>428</b>
807/014	807/016	807/018	807/023
225/014	225/016	225/018	225/023

**Pear**

<b>230</b>	<b>230B</b>	<b>231</b>	<b>232</b>	<b>366A</b>	<b>366</b>	<b>367</b>	<b>368</b>	<b>525</b>
830/010	830/012	830/014	830/016	830L/014	830L/016	830L/018	830L/021	832/033
237/010	237/012	237/014	237/016	239/014	239/016	239/018	239/021	239/033

•   •   •   •   •   •   •   •   •



**Flat End Taper (Truncated Cone)**

<b>358</b>	<b>359</b>	<b>207</b>	<b>208</b>	<b>208B</b>	<b>209</b>	<b>209B</b>	<b>365</b>	<b>235</b>	<b>236</b>	<b>237</b>	<b>210</b>	<b>210B</b>	<b>211</b>	<b>212</b>
886/016	887/018	845/009	845/010	845/012	845/014	845/016	845/021	846/012	846/014	846/016	846L/010	846L/012	846L/014	846L/016
131/016	130/018	170/009	170/010	170/012	170/014	170/016	170/021	171/012	171/014	171/016	171L/010	171L/012	171L/014	171L/016


•   •   •   •   •   •   •   •   •   •   •   •   •   •



**HP • PM • H Diamonds**



**RA • CA • W Diamonds**



Equilibrating Contouring					Wheel Round Edge			Wheel				
<b>419</b>	<b>400</b>	<b>500</b>	<b>501</b>	<b>502</b>	<b>435</b>	<b>436</b>	<b>516</b>	<b>144</b>	<b>145</b>	<b>146</b>	<b>154</b>	<b>155</b>
899/027	811/031	811/032	811/037	811/047	909/037	909/040	909/045	815/012	815/016	815/018	816/012	816/016
033/027	038/031	038/032	039/037	038/047	068/037	068/040	068/045	042/012	042/016	042/018	044/012	044/016

							Bulk Reduction				Egg		
<b>342</b>	<b>343</b>	<b>344A</b>	<b>344B</b>	<b>344C</b>	<b>344</b>	<b>446</b>	<b>511</b>	<b>512</b>	<b>513</b>	<b>514</b>	<b>255</b>	<b>396</b>	<b>416</b>
862/014	862/016	863/012	863/014	863/016	863/018	863/024	900/014	900/015	900/018	900/016	379/014	379/018	379/023
249/014	249/016	250/012	250/014	250/016	250/018	250/024	517/014	519/015	513/018	521/016	277/014	277/018	277/023

									Flame							
<b>373B</b>	<b>373C</b>	<b>373</b>	<b>374</b>	<b>408</b>	<b>373L</b>	<b>414</b>	<b>385A</b>	<b>385</b>	<b>378</b>	<b>165</b>	<b>286</b>	<b>287</b>	<b>288</b>	<b>289</b>	<b>340</b>	<b>341</b>
879K/012	879K/014	879K/016	879K/018	879K/021	879KL/016	879KL/018	885K/016	885K/016	392/016	889/009	860/010	860/012	860/014	860/016	862/010	862/012
299/012	299/014	299/016	299/018	299/021	300/016	300/018	210/016	211/016	465/016	540/009	247/010	247/012	247/014	247/016	249/010	249/012

Flame											
<b>378</b>	<b>165</b>	<b>286</b>	<b>287</b>	<b>288</b>	<b>289</b>	<b>340</b>	<b>341</b>				
392/016	889/009	860/010	860/012	860/014	860/016	862/010	862/012				
465/016	540/009	247/010	247/012	247/014	247/016	249/010	249/012				

<b>S331</b>	<b>S332</b>	<b>S333</b>	<b>S335</b>	<b>S336</b>	<b>S337</b>	<b>S341</b>	<b>S300</b>	<b>S436</b>	<b>S516</b>

JUD			
<b>241</b>	<b>296</b>	<b>300</b>	<b>530</b>
368/016	368/018	358/023	369/033
257/016	257/018	257/023	254/033



*The Right Instrument  
for Every Operation  
The High Precision  
for Unique Value*