

Dental C-bright teeth whitening bleaching light accelerator

Order No. EA1019



* Input power: 100-240V AC, 50/60Hz * Adjustable timer: 10s, 30s, 1m, 2m, 5m, 8m, 10m, 15m, 20m, 30m.

* Light resource: 6pcs high power LED

* Wavelength: 420nm-490nm * Packing dimension: 117cm*67cm*15cm * Weight: 26KG

* Light Intensity: 6 intensities, maximum >6000 mW/cm²

DENTAL SUCTION

Model No. :GS-02 - Power :1100W - Vacuum : -13KPa- Pump Flow: 600L/min - Volt./Hz:110-240V 50/60Hz



Order No. EA1020