

Amalgamator

High Speed Dental Amalgamator Dental Equipment For Mixing All Capsules With CE:

High-Tech G8 Dental Amalgamator Device Dental Equipment With CE

The machine has TÜV ISO13485, CE Certificate for suit EU standards.

HL-AH G8 Amalgamator, which controlled by a microprocessor to ensure a precise, controlled, consistent mix. With wide oscillating range (adjust any Speed as your desire from 2800 rpm to 4700rpm; time from 3S to 60S), it is suitable for all the capsule and glass ionomer sold in the market. For safe operation, the motor will stop, when the mask is open. It can satisfy different customers' requirements.



Features:

1. With digital LCD showing speed and time.
2. Adjust Time and Speed freely.
(Time from 3S~60S, Speed from 2800rpm~4700rpm)
3. Time setting with memory function.
4. Security Guard design, motor will stop when the mask is open.
5. Ergonomic design, which looks like a beetle car.
6. It's revolve according to "8" model.
7. Noiseless, very stable operation.
8. Electrical safety standard met.
9. Three colors for choose.
10. Easy maintenance.

Hand-free Impression Material Alginate Mixer for Dentistry With CE

Release your hand, Enjoy your special "car"!

HL-YMC 3 Mixer is our third generation alginate mixer, the first semi-automatic alginate mixer on the domestic market. This equipment is used for making false teeth; it replaces the handwork of mixing the Impression Material. It has characteristics of labor-saving, fast, well mix and less bubble. This new type with a high-end food control, you can adjust the speed by your foot. It will bring you excellent experience.

The unit is used to mix alginate for making dentures model.



Precision Right.



Since 1993



Main parameters:

1. Power supply: AC 220V±10%, 50 Hz±1Hz or AC 110V, 60Hz
2. Speed: you can adjust the speed as you like from 0-300rpm, just like driving car, and control your accelerator
 - A) Economical price.
 - B) Convenient operation-you can adjust the proportion of water or impression material during operation.
 - C) Four types can meet different customers' demand.
 - D) HL-YMC IV brings the story of tortoise and hare race into operation, making it vividly, and humorously.