

Air Motor Pack, 4-Hole



Features

- Package include Air Motor, Straight Head, Contra Angle and Accessories
- 20,000 rpm, Midwest type 4 hole
- Smooth Start-up, Durable, No Vibration and Quiet Operation
- Generates High Torque throughout the entire speed range
- NSK E-type Interchangeable
- It can be connected with various made attachment like KaVo, NSK, W&H and others.
- Autoclavable

Specifications

- Shank Size: 2.35 mm
- 22,000 RPM at 36 psi, 25,000 RPM at 43 psi, 27,000 RPM at 57 psi
- 42 NI/min Air Consumption at 36 psi, 51 NI/min Air Consumption at 43 psi, 72 NI/min Air Consumption at 57 psi

Contra Angle



Straight Handpiece

