

Standard 4-Hole Screw Type High-Speed Hand piece



- NSK standard screw type canister or push button canister
- 300,000-400,000 RPM
- Special silencer for quiet operation, below 66 decibels
- Fully Autoclavable
- The streamlined design feels comfortable in hand reducing fatigue
- The precision metal chuck ensures maximum concentricity of bur and proves for better wear and breakage proof, thus longer service life
- Thoroughly tested before leaving the factory

Standard 4-Hole Push Button High-Speed Hand piece



- NSK standard push button canister
- 300,000 - 400,000 RPM
- Special silencer for quiet operation, below 66 decibels.
- Fully Autoclavable
- The streamlined design feels comfortable in hand reducing fatigue
- The precision metal chuck ensures maximum concentricity of bur and proves for better wear and breakage proof, thus longer service life
- Thoroughly tested before leaving the factory

Torque 4-Hole Push Button High-Speed Handpiece



- NSK Torque push button canister
- 300,000 - 400,000 RPM
- Special silencer for quiet operation, below 66 decibels.
- Fully Autoclavable
- The streamlined design feels comfortable in hand reducing fatigue
- The precision metal chuck ensures maximum concentricity of bur and proves for better wear and breakage proof, thus longer service life
- Thoroughly tested before leaving the factory