

Curing Light



LED.B

- Wireless.
- Time setting: 5s, 10s, 15s, 20s.
- Large capacity battery. A full charge can be used for more than 500 times continuously under 10s working time mode.
- Low standby power consumption with seventy days' standby time.
- Constant light intensity. The solidification effect is not affected by the consumption of remaining power.

Power input	AC100V-240V 50Hz/60Hz
Light output	1000mW/cm ² -1200mW/cm ²
Net weight	145g
Dimensions	260mm×31mm×34mm

Centrifugal Casting Machine



Vibrator



Voltage:220V/110V
Power:70W
G.W.:2KG


Technician Polisher



Technician Polishing Machine



Hydraulic Press



Model Trimmer

