



DENTAL TECHNOLOGIES

Improving smiles around the world



ZINC PHOSPHATE

REF Code	Product
305.0000.001	Light Yellow Kit
305.0002.001	White Kit

Kit Contains:
- 32g Powder . 16.5 ml Liquid
Measuring Scoop

Order No. AZ038



Order No. AZ039

ALPHA-EASE DENSITIZER

Alpha-Ease is a Hema/Glutaraldehyde desensitizer indicated for use on sensitive cavity preparation prior to the application of restorative material. It provides a physiological seal within the dentinal tubules by polymerization with the proteins in the tubeles. The material provides virtually immediate relief from pain caused by hypersensitivity and is formulated with an antimicrobial to prevent microleakage.



Order No. AZ041

MIXING PADS

Tear-off, disposable coated sheets allow easy mixing of materials by preventing absorption into the surface. Great for mixing any two part materials including cements, impression and filling materials



Order No. AZ040

ALPHA-DENT HOMECARE FLUORIDE GEL 0.40% STANNOUS FLUORIDE

Ideal product for those patients suffering from sensitivity caused from tooth whitening procedures. The home care fluoride gel is a fluoride treatment gel intended for home use by simply brushing on and rinsing. It provides the benefits of stannous fluoride: inhibits plaque growth, helps eliminate hyper sensitivity, and aids in the prevention of dental caries. The gel is available in four great flavors! Reduces the possibility of de-calcification and band marks for orthodontic patients.



ALPHA-PRO PROPHYLAXIS PASTE

Order No. AZ042

Alpha-pro is a blend of polishing and cleaning agents containing sodium fluoride . This splatter -free formula rinses easily and completely. Medium grit is recommended for normal cleaning, polishing, and plaque removal. Coat grit is recommended for extra heavy stain, including tobacco and heavy plaque deposits. Also available without added color pigments, flavor and fluoride for patients who may allergic to additives. It is a medium grit paste ideal prior to composite bonding, sealants, and orthodontic bracket applications



ALPHA-CAINE TOPICAL ANESTHETIC

Order No. AZ043

Alpha-caine is a 20% benzocaine topical anesthetic gel used to provide fast temporary relief to patient's with minor mouth irritations, ulcers, wounds or for relieving pain prior to injection.

Four fantastic flavors:

- Strawberry
- MintBubble Gum
- Cherry