

Provis

Range of Application:

Provis is a temporary zinc oxide-based sealing compound for cavities. It has a particularly smooth and pliable consistency, optimal border sealing and a low cure time of only 15 minutes.

Adhesion is excellent even on moist dentin. This product is always easily to mold.

Provis does not irritate the dental pulp and can be easily removed with a probe or an excavator.



Order No. FP050

Steria Spray



Range of Application:

Steria Spray is a high-quality cleaning and lubricating spray handpieces, Contra angles and air motors. Steria spray creates a protective film that minimize friction and abrasion. Daily use extends the life of your sophisticated dental tools.

Order No. FP051