

**R C Help : Lubricating paste for RC Biomech**

- EDTA paste
- Contains Carbamide Peroxide
- Increases RCT instrument's life
- Gentle O2 release while irrigating with NaOCl
- Cleanses canal



**Order No. PR052**

**Available :** 3 gms Syringe with 5 Needles

**R C Cal : Intra canal dressing paste**

- Radio opaque Ca (OH)<sub>2</sub>.
- Contains barium sulphate.
- Ready to use paste.
- Highly alkaline.
- Can be easily removed by reamers.
- Non drying paste.



**Order No. PR053**

**Available :** 2 gms Syringe with 5 Needles

**R C Solve : GP & ZnOE solvent**

- Softens GP points.
- Softens zinc oxide cement.
- Controlled penetration to avoid disturbing apical seal.
- Shall not retard composite setting.



**Order No. PR054**