

## Calcium Hydroxide Powder : Dental grade

- It can be used for pulp capping. Need to mix with water for injection to make paste of desired consistency while using.

**Available** : 10 gms Bottle



## Eugenol : Dental grade

- Used with zinc oxide for temporary filling.

## Zinc Oxide : Dental grade

**Order No. PR058**

- Temporary filling material

**Available** : 100 gms Bottle



**Order No. PR059**



**Order No. PR060**

## MGS Mouthwash Tablets

- Vibrant colours • Fresh flavours • Patients choice



**Order No. PR061**

- Antiseptic effervescent tablets
- Various colours and flavours
- Freshens and cleans mouth before, during and after dental treatments
- Added detergents to effectively remove thick saliva
- Extremely economical and long shelf life