

DISPOSABLE DENTAL NEEDLES


DENJECT 덴젝트



Denject is Soft

"BIODENT CO.,LTD" are specialized to produce and export 'Dental Needles' from 1991 to now.

We acquired licence of medical device manufacture from Korea Food and Drug Administration and we have certification ISO9001,FDA and CE.



We have the high technology for 'Dental Needles'.
The DENTAL NEEDLES are designed for screw-on type syringes.
Denject's plastic hubs are pre-threaded for easy, straight attachment onto the syringe.

They are tri-bevel point, anti-coring lancet point, the round cutting edge, silicone coating and special heat treatment.
The seal of the container and the hub are colour coded to eliminate mistake or confusion.

There are 27G(35mm, 30mm, 21mm) and 30G(25mm, 21mm, 12mm) in our regular products. In addition, we have the wide range of sizes.

We have the best quality control system. Because the point is very sensitive, the control system is very important,
So we are doing our best for effective system. Also everything is worked in the clean room.



Description	Color	Size
27G XXL Long	Purple	27Gx38mm(0.41x38mm)
27G XL Long	Red	27Gx35mm(0.41x35mm)
27G Long	Orange	27Gx30mm(0.41x30mm)
27G Medium	Brown	27Gx25mm(0.41x25mm)
27G Short	Yellow	27Gx21mm(0.41x21mm)
30G Long	Violet	30Gx25mm(0.31x25mm)
30G Short	Green	30Gx21mm(0.31x21mm)
30G x.Short	Skyblue	30Gx16mm(0.31x16mm)
30G xx.Short	Blue	30Gx12mm(0.31x12mm)

ISO9001 0499