



**Order No. MA022**

### Alginplus

Extra-high precision alginate for impressions with chromatic phase indicator.



**Order No. MA023**

### Ormament Putty

Polysiloxane based, condensation-curing material for preliminary precision impressions.



**Order No. MA024**

### Ormamax Light

Polysiloxane based, condensation-curing material for second precision impressions.



**Order No. MA025**

### Ormamax Regular

Polysiloxane based, medium viscosity condensation-curing material for precision impressions.



**Order No. MA026**

### Ormactivator Gel

Gel catalysts for condensation-curing silicones.